

Microbiology and Coagulation Testing: LIS


Pathology Informatics Summit 2018

May 21, 2018

Dr. JiYeon Kim

Goals of talk

1. Review what makes for distinctive microbiology LIS features, and what data drives decision-making in clinical infectious diseases
2. Be familiar with trends in microbiology automation, new diagnostic technologies
3. Understand how new anticoagulation drugs are affecting coagulation test results
4. Consider future LIS data needs, including the need for “talking” interfaces with other clinical information systems outside the lab



Provider questions for clinical microbiology lab:


1. Is my patient's illness caused by a microbe?
2. If yes, what is it?
3. What is the antibiotic susceptibility profile of that organism so that therapy can be targeted?

Uncertainty drives use of broad-spectrum antibiotics.



Role of Microbiologist: Cartoon by Czichos

Global rise in antibiotic consumption



THE RESISTANCE MOVEMENT

Carbapenem-resistant Enterobacteriaceae have been on the move since at least 1996.



1 2000: Analysis of a 1996 sample from a North Carolinian hospital finds infectious *Klebsiella pneumoniae* carrying a gene called KPC that confers resistance to carbapenems.

2 2003: KPC-positive bacteria are found spreading rapidly through hospitals across New York City. By 2007, 21% of *Klebsiella* in the city carry the resistance gene.

3 2005: KPC-positive bacteria make their way from New York to several other countries, including Israel. From Israel, the bacteria travel to Italy, Colombia, the United Kingdom and Sweden.

4 2008: Doctors in Sweden find a new carbapenem-resistance gene, NDM. Traced back to India, NDM-positive bacteria have moved quickly.

Mortality from Infectious Diseases


Deaths due to **infectious diseases** in the United States:


1. Heart disease: 633,842
2. Cancer: 595,930
3. **Chronic lower respiratory diseases: 155,041**
4. Accidents (unintentional injuries): 146,571
5. Stroke (cerebrovascular diseases): 140,323
6. Alzheimer's disease: 110,561
7. Diabetes: 79,535
8. **Influenza and Pneumonia: 57,062**
9. Nephritis, nephrotic syndrome and nephrosis: 49,959
10. Intentional self-harm (suicide): 44,19
11. **Septicemia: 40,685**

DRIVING INCREASES:

- Aging population
- Drug-resistance microorganisms
- Global warming, spread of tropical disease carriers into new areas

Global Clinical Microbiology Market







Microbiology LIS

Interfaces



Role of Microbiologist: Cartoon by Czichos






Microbiology LIS

Interfaces

- 2 Programmable plate inoculators/streakers
- 3 Automated incubators; automated growth and detection systems
 - Ex: Bactec, BacT/ALERT, MGIT, TREK, etc.
- 4 Automated plate readers

GeneXpert




Xpert SA Nasal Complete G3

Assay 4

Test Result: **MRSA POSITIVE; SA POSITIVE**

Test and Analyte Result

Analyte Name	Ct	EndPt	Analyte Result	Probe Check Result
SPC	36.9	23.0	NA	PASS
SPA	23.9	428.0	POS	PASS
mec	24.4	349.0	POS	PASS
SCC	25.5	440.0	POS	PASS



Microbiology LIS

Interfaces


- 2 Programmable plate inoculators/streakers
- 3 Automated incubators; automated growth and detection systems
 - Ex: Bactec, BacT/ALERT, MGIT, TREK, etc.
- 4 Automated plate readers
- 5 5 Identification of organisms
 - 5 Biochemical-based microbial identification
 - 5 MALDI-ToF mass spectrometry
 - 5 Genetic markers and sequencing
- 6 6 Susceptibility testing
 - 6 Traditional phenotypic tests
 - 6 Genetic markers

Specimen types, and related preparation and sequences for testing, are varied and difficult to completely automate. Results take time (hours to days).


Microbiology LIS

Best Practices

- ✓ Intelligent decision support business rules for routine common tasks (efficiency) as well as sentinel events
- ✓ Preliminary versus final status messages.



Automation and New Tech



Data is relational and needs traceability:

Order (e.g., respiratory culture)

→ multiple **specimens** (e.g., induced sputum #1, #2, #3)


→→ multiple **culture media** (e.g., aerobic, anaerobic, fungal)

→→→ multiple **isolates** (e.g., Staph. aureus, H. influenzae)

→→→→ multiple **observations** per isolate (e.g., gram stain, colony count, identification, antibiotic **susceptibility** panels), separated as “preliminary” versus “final”

Microbiology LIS

Best Practices



Generic Microbiology HL7 Result Message		Example of sputum culture
MSH	Message header	
PID	Patient identification	
[PV1]	Patient visit information	
{	--- Order begins	
[ORC]	Order control	
OBR	Observation request	<u>OBR 1 = sputum culture</u> OBX 1 = Microorganism 1 = S. aureus OBX 2 = Colony ct 1 = 10,000-90,000 OBX 3 = Microorganism 2 = Beta-hemolytic Streptococcus OBX 4 = Colony ct 2 = <1,000 OBX 5 = Microorganism 3 = Haemophilus influenzae OBX 6 = Colony ct 3 = 10,000-90,000
[{ NTE }]	Notes and comments	
[OBX]	Observation <i>OBX-4-observation sub-ID sequential numbering to link observations to isolates</i>	
[{ NTE }]	Notes and comments	
}	--- Observation end	
}	--- Order end	

Generic Microbiology HL7 Result Message		Example of sputum culture
MSH	Message header	<u>OBR 1 = sputum culture</u>
PID	Patient identification	OBX 1 = Microorganism 1 = S. aureus
[PV1]	Patient visit information	OBX 2 = Colony ct 1 = 10,000-90,000
{	--- Order begins	OBX 3 = Microorganism 2 = Beta-hemolytic Streptococcus
[ORC]	Order control	OBX 4 = Colony ct 2 = <1,000
OBR	Observation request	OBX 5 = Microorganism 3 = Haemophilus influenzae
[{ NTE }]	Notes and comments	OBX 6 = Colony ct 3 = 10,000-90,000
[OBX]	Observation <i>OBX-4-observation sub-ID sequential numbering to link observations to isolates</i>	<u>OBR 2 = susceptibility panel for S aureus</u>
		OBX 1 = Ampicillin MIC 1 = 32 ug/mL = R
		OBX 2 = Amox+Clav MIC 1 = 2 ug/mL = S
		OBX 3 = Cefazolin MIC 1 = 8 ug/mL = S
		OBR 4 = Ampicillin MIC 1 = 32 ug/mL = R
[{ NTE }]	Notes and comments	<u>OBR 3 = susceptibility panel for H influenzae</u>
}	--- Observation end	OBX 1 = Ampicillin Kirby-Bauer 3 = 3 = S
}	--- Order end	OBX 2 = Amox+Clav Kirby-Bauer 3 = 3 = S
		OBX 3 = Cefazolin Kirby-Bauer 3 = 3 = S


```

MSH|^~\&|LabOneApp|LabOne^45D0470381^CLIA|NEDSS^1644^WA-
DOH||20011001183345||ORU^R01^ORU_R01|113661|P|2.3.1<cr>
PID|||999-3^^^GoodDr^MR||Able^Michael^D^^^^L<cr>
ORC|RE|0889436^GoodDr|ABC012345^LabOne|||20011001182914|
L0234^Roberts^Steve^^^^^LabOne^L^^^EI
|^Good^Robert^^^^MD^^L<cr>
OBR|1|0889436^GoodDr|ABC012345^LabOne|6460-0^Spt Routine
Cult^LN|||20011001091234|||200110010823|SPT&
Sputum&HL70070|^Good^Robert^^^^MD^^L|||20011002072359
||MB|P<cr>
OBX|1|CE|11475-1^MICROORGANISM IDENTIFIED:^LN|1|L-
24801^Staphylococcus aureus^SNM<cr>
OBX|2|CE|564-5^Colony count^LN|1|10,000-90,000<cr>
OBX|3|CE|11475-1^MICROORGANISM IDENTIFIED:^LN|2|L-
25128^Beta hemolytic Streptococcus A^SNM<cr>
OBX|4|CE|564-5^Colony count^LN|2|<1,000<cr>
OBX|5|CE|11475-1^MICROORGANISM IDENTIFIED:^LN|3|L-13401 ^
Haemophilus influenzae^SNM<cr>
OBX|6|CE|564-5^Colony count^LN|3|10,000-90,000<cr>
OBR|2||ABC012346^LabOne|29576-6^Bacterial Susc Panel
IsIt^LN||||||||||||||||11475-1&MICROORGANISM
IDENTIFIED:&LN^1^Staphylococcus
aureus||0889436&GoodDr^ABC012345&LabOne<cr>
OBX|1|CE|28-1^Ampicillin MIC^LN|1|32|µg/mL| |R<cr>
OBX|2|CE|32-3^Amoxicillin+Clav MIC^LN|1|2|µg/mL| |S<cr>
OBX|3|CE|76-0^Cefazolin MIC^LN|1|8|µg/mL| |S<cr>
OBR|3||ABC012347^LabOne|29576-6^Bacterial Susc Panel
IsIt^LN||||||||||||||||11475-1&MICROORGANISM
IDENTIFIED:&LN^3^Haemophilus
influenzae||0889436&GoodDr^ABC012345&LabOne<cr>
OBX|1|CE|29-9^Ampicillin KB^LN|3| |||S<cr>
OBX|2|CE|21-6^Amoxicillin+Clav KB^LN|3| |||S<cr>
OBX|3|CE|77-8^Cefazolin KB^LN|3| |||S<cr>

```

Unlike chemistry tests (close-ended, from 1 order to 1 result), microbiology data requires flexibility to accommodate many **open-ended parent-child relationships** (1 to many).

Specimen types, and related preparation and sequences for testing, are varied and difficult to completely automate. Results take time (hours to days)

Data is relational and n

Order (e.g., respiratory cult
→ multiple **specimens** (e.g.,
→→ multiple **culture** medi
→→→ multiple **isolates** (e
→→→→ multiple **observa**
identification, antibiotic sus
“final”

! URINE CULTURE

Status: Final result Visible to patient: No (Not Released) Next appt: None Dx:

Order:

Notes Recorded by

4/11/18 12:00 PM

Flag A

FINAL RESULT

! (A)
>100,000 cfu/ml Escherichia coli
Organism is a possible multi-drug resistant bacteria.
In serious infections cefepime may not be as effective as other agents. Please consult Infectious Diseases regarding its use.
<10,000 CFU/mL of Insignificant Growth

ORGANISM

Escherichia coli !

A

Susceptibility

	Escherichia coli	
	ANTIMICROBIAL SUSCEPTIBILITY	
Ampicillin	>=32	RESISTANT
Cefazolin urine	>=64	RESISTANT
Cefepime	4	INTERMEDIATE
Ceftazidime	16	RESISTANT
Ceftriaxone	>=64	RESISTANT
Ciprofloxacin	>=4	RESISTANT
Ertapenem	<=0.5	SUSCEPTIBLE
Gentamicin	<=1	SUSCEPTIBLE
Meropenem	<=0.25	SUSCEPTIBLE
Nitrofurantoin	<=16	SUSCEPTIBLE
Piperacillin+Tazobactam	<=4	SUSCEPTIBLE
Trimethoprim+Sulfamethoxazole	<=20	SUSCEPTIBLE

Narrative

Specimen Collected: 04/11/18 12:00 PM

Last Resulted: 04/14/18 2:22 AM

Lab Flowsheet Order Details View Encounter Lab and Collection Details Routing Result History

ctices

business rules for routine
well as sentinel events

s messages.


, OBX-4-observation sub-ID
observations to isolates

g sent to EHR (if it can

Outbreaks: Hospitals, Communities

Data needs aggregation, mandatory reports to public health authorities

- from **multiple specimens** to **one patient** (e.g., same pathogen isolated from multiple sources)
- from **multiple antibiotic susceptibility patterns** from the same and/or different organisms to the **appropriate treatment plan**
- from **multiple patients** to the **same hospital floor or unit**
- ... to the **same institution/city/state** (e.g., antibiogram)



Specimen types, and related preparation and sequences for testing, are varied and automate. Results take time

Data is relational:

Order (e.g., respiratory culture)
→ multiple **specimens** (e.g., ind)
→→ multiple **culture** media (e.g.,)
→→→ multiple **isolates** (e.g.,)
→→→→ multiple **observation**
identification, antibiotic suscep
“final”

Data needs aggregation, r health authorities

- from **multiple specimens** to c
from multiple sources)
- from **multiple antibiotic suscep**
different organisms to the ap
- from **multiple patients** to the
- ... to the **same institution/city**

! AEROBIC BACTERIAL CULTURE, STERILE SITE

Status: Final result Visible to patient: No (Not Released) Next appt: 05/18/2018 at 10:20 AM in Internal Medicine

Newer results are available. Click to view them now.

FINAL RESULT	8/12/17 10:46 AM	Flag A
! (A) Moderate growth of: Enterobacter cloacae Complex Moderate growth of: Klebsiella oxytoca Organism is a possible multi-drug resistant bacteria. In serious infections cefepime may not be as effective as other agents. Please consult Infectious Diseases regarding its use. Moderate growth of: Enterococcus faecalis		
ORGANISM	Enterobacter cloacae complex !	A
ORGANISM	Klebsiella oxytoca !	A
ORGANISM	Enterococcus faecalis !	A
Susceptibility		
	Enterobacter cloacae complex ANTIMICROBIAL SUSCEPTIBILITY	Klebsiella oxytoca ANTIMICROBIAL SUSCEPTIBILITY
Ampicillin	>=32 RESISTANT	>=32 RESISTANT
Ampicillin+Sulbactam		>=32 RESISTANT
Cefazolin	>=64 RESISTANT	>=64 RESISTANT
Cefepime		2 SUSCEPTIBLE
Ceftazidime		16 RESISTANT
Ceftriaxone	<=1 SUSCEPTIBLE	>=64 RESISTANT
Ciprofloxacin	<=0.25 SUSCEPTIBLE	<=0.25 SUSCEPTIBLE
Ertapenem		<=0.5 SUSCEPTIBLE
Gentamicin	<=1 SUSCEPTIBLE	<=1 SUSCEPTIBLE
Meropenem		<=0.25 SUSCEPTIBLE
Piperacillin+Tazobactam	<=4 SUSCEPTIBLE	8 SUSCEPTIBLE
Trimethoprim+Sulfamethoxazole	<=20 SUSCEPTIBLE	>=320 RESISTANT
Susceptibility		
	Enterococcus faecalis ANTIMICROBIAL SUSCEPTIBILITY	
Ampicillin	<=2 SUSCEPTIBLE	
Gentamicin Synergy	SUSCEPTIBLE ¹	
Streptomycin Synergy	SUSCEPTIBLE ¹	
Vancomycin	1 SUSCEPTIBLE	

Practices

Business rules for routine
as well as sentinel events
status messages.



ure, OBX-4-observation sub-ID
link observations to isolates

being sent to EHR (if it can
b”

reports

lic health

New species, larger microbial world



Nature Reviews | Microbiology

Specimen types, and related preparation and sequences for testing, are varied and difficult to completely automate. Results take time (hours to days).

Data is relational:

Order (e.g., respiratory culture)

→ multiple **specimens** (e.g., induced sputum #1, #2, #3)

→→ multiple **culture** media (e.g., aerobic, anaerobic, fungal)

→→→ multiple **isolates** (e.g., Staph. aureus, H. influenzae)

→→→→ multiple **observations** per isolate (e.g., gram stain, colony count, identification, antibiotic susceptibility panels), separated as “preliminary” versus “final”

Data needs aggregation, mandatory reports to public health authorities

- from multiple specimens to one patient (e.g., same pathogen isolated from multiple sources)
- from multiple antibiotic susceptibility patterns from the same and/or different organisms to the appropriate treatment plan
- from multiple patients to the same hospital floor or unit
- ... to the same institution/city/state (e.g., antibiogram)

Results and result categories evolve

- **taxonomy** frequently gets updated (names change) – multiple historic names for equivalent organisms
- **susceptibility breakpoints** change (what is defined as susceptible or resistant)

Microbiology LIS

Best Practices


- ✓ Intelligent decision support business rules for routine common tasks (efficiency) as well as sentinel events
- ✓ Preliminary versus final status messages.
- ✓ HL7 2.5.1 OBR/OBX structure, OBX-4-observation sub-ID sequential numbering to link observations to isolates
- ✓ Support structured data being sent to EHR (if it can receive it) versus “text blob”
- ✓ Support viewing related reports
- ✓ Required reporting to public health
- ✓ Update names frequently and connect historic to current names, manage changing interpretation criteria
- ✓ Use of standard nomenclature with LOINC, SNOMED

Faster, smarter

Provider questions for clinical microbiology lab:

1. Is my patient's illness caused by a microbe?
2. If yes, what is it?
3. What is the antibiotic susceptibility profile of that organism so that therapy can be targeted?

Uncertainty drives use of broad-spectrum antibiotics.





3:53

+ QUEUE

DOWNLOAD

EMBED

TRANSCRIPT



PUBLIC HEALTH


Federal Efforts To Control Rare And Deadly Bacteria Working

April 3, 2018 - 3:29 PM ET

Heard on *All Things Considered*



RICHARD HARRIS





The CDC is trying to stop E. coli and other bacteria that have become resistant to antibiotics because they can cause a deadly infection.

Science Photo Library/Getty Images/Science Photo Library

Testing... included **carbapenemase production testing** and **molecular detection of genes** encoding for the five carbapenemases of primary public health concern... laboratories were asked to report positive findings to local public health authorities and CDC **within 1 day**...

For each carbapenemase-producing isolate detected... state health department staff members contact the health care facility to **review infection control measures** and consider performing **on-site infection control assessments**. If indicated, **contacts of the index patient are screened** to detect transmission... Response activities continue until transmission is **controlled**.

- *The percentage of ESBL phenotype Enterobacteriaceae decreased by 2% per year; by comparison, the CRE percentage decreased by 15% per year.*





3:53

+ QUEUE

DOWNLOAD

EMBED

TRANSCRIPT



PUBLIC HEALTH


Federal Efforts To Control Rare And Deadly Bacteria Working

April 3, 2018 - 3:29 PM ET

Heard on *All Things Considered*



RICHARD HARRIS





The CDC is trying to stop E. coli and other bacteria that have become resistant to antibiotics because they can cause a deadly infection.

Science Photo Library/Getty Images/Science Photo Libra

• The proportion of *Enterobacteriaceae* infections that were CRE remained lower and decreased more over time than the proportion that were ESBL phenotype. This difference might be explained by the **more directed control efforts implemented** to slow transmission of CRE than those applied for ESBL-producing strains.

• These data suggest that an **early aggressive response**, as outlined in CRE-specific infection prevention recommendations released beginning in 2009, can **slow emergence and even decrease the occurrence of infections** from resistant pathogens.



Faster, smarter



Xpert SA Nasal Complete G3				
Assay 4				
Test Result: MRSA POSITIVE; SA POSITIVE				
Test and Analyte Result				
Analyte Name	Ct	EndPt	Analyte Result	Probe Check Result
SPC	36.9	23.0	NA	PASS
SPA	23.9	428.0	POS	PASS
mec	24.4	349.0	POS	PASS
SCC	25.5	440.0	POS	PASS

Removing patients from unnecessary contact precautions with single negative PCR compared to 3 consecutively negative MRSA cultures from nasal swabs


Strategy	Passive cultures	Active surveillance	PCR screening
Discontinuation rates of contact precautions	6.6%	26.2%	63.8%
Fewer contact precaution days	104	418	1,841
Cost savings	\$86,950	\$349,472	\$1,539,180

Faster, smarter

Provider questions for clinical microbiology lab:

1. Is my patient's illness caused by a microbe?
2. If yes, what is it?
3. What is the antibiotic susceptibility profile of that organism so that therapy can be targeted?

We have actionable data. We need to develop smart message interfaces from the LIS to the **pharmacy ...
From the LIS to **ADT** for bed management...**



Faster, smarter


Provider questions for clinical microbiology lab:

1. Is my patient's illness caused by a microbe?
2. If yes, what is it?
3. What is the antibiotic susceptibility profile of that organism so that therapy can be targeted?

Is it more important to identify organisms or identify resistance quickly?

How to clarify colonization from infection? (Qty?)


How to ingest and report microbiome data? Which are the “good” versus “bad”? Is it important to report both?



Coagulation LIS


Prothrombin Time

Measures time of fibrin clot formation, beginning with factor VII activation, through the final step in which fibrinogen is converted to fibrin



Common Pathway

Extrinsic Pathway – PT



Activated Partial Thromboplastin Time

Intrinsic Pathway – aPTT

Activator (e.g., silica, kaolin, celite, elegaic acid)

XII

XIIa

XI

XIa

IX

Ca²⁺

IXa
VIIIa

X

Ca²⁺
PL

Xa
Va

IIa

(thrombin)

II

(prothrombin)


Common Pathway

fibrin
clot

I

(fibrinogen)

Measures time it takes to form a fibrin clot, beginning with factor XII activation, through the final step in which fibrinogen is converted to fibrin



Clot-based tests

Intrinsic Pathway – aPTT

- aPTT
- Protein C
- Protein S
- Lupus anticoagulant
- Factor VIII
- Factor IX
- Factor XI
- Factor XII

HMWK, Kallikrein

XII

PK → Kallikrein

XIIa

HMWK

XI

XIa

Ca²⁺

IX

IXa

VIIIa

Ca²⁺

Extrinsic Pathway – PT

- PT
- INR
- Factor II
- Factor V
- Factor VII
- Factor X

TF (in PL)

Ca²⁺

VIIa/TF

VII

PL, Ca²⁺

Xa

Va

PL, Ca²⁺

IIa

(thrombin)

II

(prothrombin)

V, VIII

D-dimer
FDPs
Common Pathway


fibrin
clot

(fibrinogen)

Thrombomodulin
(Endothelial cells)

Protein C,
Protein S

- Thrombin time



New oral anticoagulants

Initial U.S. Approval: 2012


Eliquis
(apixaban) tablets 500mg

This site is intended for U.S. residents 18 years of age or older.

Call Us at 1.855.ELIQUIS

Activate Your Co-pay Card


Tell Me How ELIQUIS Can Help



Reduce the Risk of Stroke Due to Nonvalvular Atrial Fibrillation (AFib)

ELIQUIS & Nonvalvular AFib


For people with AFib, a type of irregular heartbeat, not caused by a heart valve problem



Treat Blood Clots In the Veins of the Legs/Lungs & Reduce Risk of it Occurring Again

ELIQUIS & DVT/PE

For people with deep vein thrombosis (DVT) or pulmonary embolism (PE)



Decrease Risk of DVT Blood Clots After Hip or Knee Replacement Surgery

ELIQUIS & Hip/Knee Replacement

For people who have had hip or knee replacement surgery

SELECTED IMPORTANT SAFETY INFORMATION & INDICATIONS

For people taking ELIQUIS® (apixaban) for atrial fibrillation: Do not stop taking ELIQUIS without talking to your doctor. Stopping ELIQUIS increases your risk of having a stroke.

Pradaxa
dabigatran etexilate
CAPSULES

For U.S. Residents Only

FOR HEALTHCARE PROFESSIONALS

Important Safety Information

Prescribing Information

Medication Guide

Call for support 1-877-481-5332

Treating with a Blood Thinner

What is PRADAXA?

Starting PRADAXA

Staying on PRADAXA

Caregivers

Resources

SIGN UP for FREE Support

If you have atrial fibrillation (AFib) not caused by a heart valve problem, PRADAXA is proven to lower your risk of stroke and is the only blood thinner, other than warfarin, with a **specific reversal treatment**. Available in over 3,100 hospitals nationwide.*

SEE how PRADAXA

DOWNLOAD the

IMPORTANT SAFETY INFORMATION AND USE OF PRADAXA

SEE MORE

Do not stop taking PRADAXA without talking to the doctor who prescribes it for you. PRADAXA may need to be stopped prior to surgery or a medical or dental procedure. When you stop taking PRADAXA, your blood may clot from forming.

Initial U.S. Approval: 2014

Xarelto
rivaroxaban tablets

Home About XARELTO® Health Conditions Programs and Stories Cost Support and More

XARELTO® or aspirin? There's more to know about reducing the risk of another DVT or PE. [Learn More](#)

Once they got the facts, they chose XARELTO®, a latest-generation blood thinner.

Learn how XARELTO® works

Janssen CarePath

IMPORTANT SAFETY INFORMATION

WHAT IS THE MOST IMPORTANT INFORMATION I SHOULD KNOW ABOUT XARELTO®?

For people taking XARELTO® for atrial fibrillation: People with atrial fibrillation (an irregular heart beat) are at

WHAT IS XARELTO®?

XARELTO® is a prescription medicine used to reduce the risk of stroke and blood clots in people with atrial fibrillation, not caused by a heart valve problem. For

Initial U.S. Approval: 2017

Clot-based tests

Intrinsic Pathway – aPTT

- aPTT
- Protein C
- Protein S
- Lupus anticoagulant
- Factor VIII
- Factor IX
- Factor XI
- Factor XII

HMWK, Kallikrein

XII

PK → Kallikrein

XIIa

HMWK

XI

XIa

Ca²⁺

IX

IXa

VIIIa

Extrinsic Pathway – PT

- PT
- INR
- Factor II
- Factor V
- Factor VII
- Factor X

TF (in PL)

Ca²⁺

VIIa/TF

Ca²⁺

VII

Ca²⁺

PL, Ca²⁺

Xa

Va

PL, Ca²⁺

IIa

(thrombin)

II

(prothrombin)

V, VIII

D-dimer
FDPs
Common Pathway

fibrin
clot

(fibrinogen)

Thrombomodulin
(Endothelial cells)

Protein C,
Protein S


Xa inhibitor:

- Rivaroxaban
- Apixaban

IIa inhibitor:

- Dabigatran


- Thrombin time



Clot-based tests: “open-ended” cascades

Intrinsic Pathway – aPTT

Extrinsic Pathway – PT



Effect of Various Anticoagulants On Commonly Used Coagulation Assays.


Source: Dorothy Adcock MD, Chief Medical Officer, LabCorp Diagnostics

Coagulation Assays	VKA (Influence)	UFH (Influence)	LMWH (Influence)	Rivaroxaban (Influence)	Apixaban (Influence)	Edoxaban (Influence)	Dabigatran (Influence)
PT	↑	no effect or ↑ ¹	=	↑	no effect or ↑ _(weak)	↑ _(weak)	↑
aPTT	↑	↑	no effect or ↑ _(weak)	↑	no effect or ↑ _(weak)	↑ _(weak)	↑
Fibrinogen (Clauss Method)	=	↓	=	=	=	=	= / ↓
Thrombin Time TT	=	↑	↑	=	=	=	↑
Factor Assays (clotting assays)	↓	aPTT based: ↓ ² PT based: =	aPTT based: ↓ ² PT based: =	↓ ²	↓ ² / = _(weak)	↓ ²	↓ ²
DDi, VWF: Ag, VWF: RCo	=	=	=	=	=	=	=
Anti-Xa Activity (UFH or LMWH)	=	=	=	=	=	=	=
Antithrombin Activity FXa-based Assay	=	=	=	=	=	=	=
Antithrombin Activity FIIa-based Assay	=	=	=	=	=	=	=
Protein C Activity Clot-based Assay	=	=	=	=	=	=	=
Protein C Activity Chromogenic Assay	=	=	=	=	=	=	=
Protein S Activity Clot-based Assay	=	=	=	=	=	=	=
Free Protein S Ag (Immunological Assay)	↓	=	=	=	=	=	=
Lupus Anticoagulant Testing: "sensitive" aPTT and dRVVT (screening, mixing, confirmation)	↑ ⁴	↑ ⁴	=	↑ ⁴	↑ ⁴	↑ ⁴	↑ ⁴
Resistance to Activated Protein C	↑ ⁵	↑ ¹	=	↑ ³	↑ ³	↑ ³	↑ ³
Reptilase Time	=	=	=	=	=	=	=

40% false-positive lupus anticoagulant test results at KPSC – new “indeterminate” result category; using chart review to check medication history for some cases

LIS needs Pharmacy and Vice Versa

- Many potential areas for connecting and improving patient care:
 - Oncology and pharmacogenetic testing
 - Antibiotic stewardship activities
 - Therapeutic drug monitoring (e.g., methotrexate)
 - Opioid / pain management drug prescription monitoring (adherence)
 - Safe renal dosing of medications
 - Anticoagulation interference in diagnostic work-up of hypercoagulability conditions
 - Etc...
- Some pharmacy systems are able to receive laboratory results, enabling rules-based review (e.g., therapeutic drug monitoring)
- However, in most cases, pharmacy and laboratory systems rarely communicate



<http://yuanyuanliang.blogspot.com>

Thank you!

References

- Aller, RD and Salazar, V. "Microbiology Laboratory Information Systems." Manual of Commercial Methods in Clinical Microbiology, International Edition, Second Edition. Ed. Allan L. Truant. John Wiley & Sons, Inc. 2016.
- Almonacid DE, Kraal L, Ossandon FJ, Budovskaya YV, Cardenas JP, Bik EM, Goddard AD, Richman J, Apte ZS. 16S rRNA gene sequencing and healthy reference ranges for 28 clinically relevant microbial taxa from the human gut microbiome. PLoS One. 2017 May 3;12(5):e0176555.
- CDC. Implementation Guide for Transmission of Microbiology Result Messages as Public Health Information Using Version 2.3.1 of the Health Level Seven (HL7) Standard Protocol. May 27, 2003.
- Croxatto A, Dijkstra K, Prod'hom G, and Greub G. Comparison of Inoculation with the Inoqula and WASP Automated Systems with Manual Inoculation. J Clin Microbiol. 2015 Jul; 53(7): 2298–2307.
- Fournier PE, Drancourt M, Colson P, Rolain JM, La Scola B, Raoult D. Modern clinical microbiology: new challenges and solutions. Nat Rev Microbiol. 2013 Aug;11(8):574-85.
- Laxminarayan R, et al. Antibiotic resistance-the need for global solutions. Lancet Infect Dis. 2013 Dec;13(12):1057-98.
- Rhoads DD, Novak SM, Pantanowitz L. A review of the current state of digital plate reading of cultures in clinical microbiology. J Pathol Inform 2015;6:23.
- Shenoy ES, Kim JY, Rosenberg ES, Cotter JA, Lee H, Walensky RP, Hooper DC. Discontinuation of Contact Precautions for Methicillin-Resistant Staphylococcus aureus: A Randomized Controlled Trial Comparing Passive and Active Screening With Culture and Polymerase Chain Reaction. CID 2013 Jul;57(2):176-84
- Schiff GD, Klass D, Peterson J, Shah G, Bates DW. Linking laboratory and pharmacy: opportunities for reducing errors and improving care. Arch Intern Med. 2003 Apr 28;163(8):893-900.
- Woodworth KR, Walters MS, Weiner LM, et al. Vital Signs: Containment of Novel Multidrug-Resistant Organisms and Resistance Mechanisms — United States, 2006–2017. MMWR Morb Mortal Wkly Rep 2018;67:396-401.

Clot-based tests: “open-ended”

Intrinsic Pathway – aPTT

Extrinsic Pathway – PT

